#### **AF1 SERIES**

### SINGLE CHANNEL AUTONOMOUS PUMP

ELVEFLOW.COM/MICROFLUIDIC-FLOW-CONTROL-PRODUCTS/FLOW-CONTROL-SYSTEM/HIGH-ACCURACY-PRESSURE-PUMPS/



# AN AUTONOMOUS PUMP DESIGNED TO MATCH THE NEEDS OF ALL MOBILE SCIENTISTS



The AF1 is a high performance autonomous pressure and flow controller. It comes in three different ranges and embeds pressure and vacuum sources. It is compatible with Elveflow software.

| / | STA | NDA | INNE  | UNIT |
|---|-----|-----|-------|------|
|   | JIA | NUA | LUIIL | UITI |

**✓ NO COMPUTER NEEDED** 

#### **UNIQUE PERFORMANCES**


- > Pressure resolution 100 μbar
- > Pressure stability 100 μbar
- > Response time **50 ms**
- > Settling time down to 100 ms

#### **APPLICATIONS**

- > Digital microfluidics: micro-droplets, anisotropic articles, double emulsions generation & handling
- > Beads and particles manipulation
- > Fast liquid sample switching
- > Cell culture experiments under medium perfusion

HOW IT WORKS

AF1 SERIES



## 1 AF1 Pressure and vacuum generator

Fast and accurate pressure and vacuum control for your system.

2 Sample

Push and pull liquid from the reservoir to your microfluidic chip

3 Microfluidic device

The liquid is smoothly and precisely injected into the microfluidic

4 Monitoring

Control pressure and vacuum using your computer or the instrument's front panel dial.

#### **FEATURES & BENEFITS**



Short settling time

Piezo technology allowing a blazing fast flow change in any microdevice High flow stability

Pressure stability down to 100 µbar ensuring a superior flow performance over a large flow range

Accurate flow control

Input a flow value into the software. Flow regulation down to 7.5 nL/min



Software automation

Control all instruments through a single panel. Powerful script module to automate control and injection over days Create your own program •

Software Development Kits (C++, Python, MATLAB® and LabVIEW® libraries)

Enhanced data saving

Up to 10 ms sampling rate to take out the best of your results



Easy to install and use

Start out of the box and set everything up within minutes

Several pressure range

Choose among the three pressure ranges available

Knob pressure control

Monitor and control pressure using the front panel knob and screen

| AF1 PREMIUM                       | AF1 200   | AF1 1,600                                 | AF1 DUAL  |  |
|-----------------------------------|---|---|---|--|
| Pressure range                    | <b>0 to 200 mbar</b><br>(0 to 2.9 psi)  | <b>0 to 1,600 mbar</b><br>(0 to 23 psi)   | <b>-700 to 1,000 mbar</b> (-10 to 14 psi)                                 |  |
| Type of pressure                  | positive  | positive                                  | negative & positive   |  |
| Pressure sensor resolution        | <b>0.006 % FS</b><br>12.2 μbar (0.0007 psi)   | <b>0.006 % FS</b><br>122 μbar (0.007 psi) | <b>0.006 % FS</b><br>122 µbar (0.007 psi)                                 |  |
| Pressure stability <sup>(1)</sup> | <b>100 μbar</b><br>0.05 % FS (0.0014 psi)   | <b>1 mbar</b><br>0.05 % FS (0.014 psi)    | <b>0.006 % FS</b><br>122 µbar (0.007 psi)                                 |  |
| Response time <sup>(2)</sup>      |   | 50 ms                                     |   |  |
| Settling time (3)                 | down to 100 ms  |   |   |  |
| Supply pressure<br>(min - max)    | integrated pump<br>no pressure source needed  |   | integrated pressure & vacuum source<br>no pressure & vacuum source needed |  |
| Liquid compatibility              | any aqueous or organic solvent, oil, or biological sample solution can be propelled |   |   |  |
| Output connectors                 | stainless steel female luer lock  |   |   |  |

Non-contractual information, may be changed without notice.

POWER CONSUMPTION: 15 W (100 V to 240 V - 50 Hz to 60 Hz) CASE DIMENSIONS (length x width x height): 220 x 130 x 130 mm WEIGHT: 1.7 kg TTL TRIGGER: input/output 5 V

(1) Output stability measured at 150 mBar with an external High accuracy pressure sensor (Druck DPI150) (2) Depending on user computer operating system (3) Volume dependent – Measurement done on 12 mL tank for a set point from 0 to 200 mbar

#### **PRODUCTS & SERVICES**

| ELEMENTS PROVIDED BY ELVEFLOW  | INCLUDED | OPTIONAL |
|--|----------|----------|
| Software & libraries Controll all Elveflow Instruments with the same smart interface   | •        |          |
| AF1 connection kit A complete set of accessories fitted for the AF1 pressure generator |          | •        |
| Reservoirs Gas tight reservoirs with ergonomic fluidic connection                      |          | •        |
| Flow sensors A line of sensors to monitor very low liquid flow rates                   |          | •        |
| Service The Elveflow expertise & support to offer you individually tailored solutions  | •        |          |

Non-contractual information, may be changed without notice